

# Radiation Protection Gloves

Our Radiation Protection Gloves provide superior tactile feedback without compromising attenuation.



## Barrier Radiation Protection Gloves

LEAD FREE | POWDER FREE

### XR1

|        |        |         |
|--------|--------|---------|
| 60 KVP | 80 KVP | 100 KVP |
| 46%    | 35%    | 30%     |

- .20" thick textured finger/palm
- High tensile strength
- Direct exposure attenuation rates

### XR2

|        |        |         |
|--------|--------|---------|
| 60 KVP | 80 KVP | 100 KVP |
| 64%    | 54%    | 48%     |

- .35" thick textured finger/palm
- Maximum attenuation
- Direct exposure attenuation rates



## WRP Radiaxon® Gloves

LEAD FREE | POWDER FREE

|        |        |         |
|--------|--------|---------|
| 60 KVP | 80 KVP | 100 KVP |
| 59%    | 49%    | 42%     |

- .35" thick textured finger/palm
- High attenuation and tactile sensitivity
- Direct exposure attenuation rates



## IB Radiation Protection Gloves

LEAD FREE | POWDER FREE | LATEX FREE

|        |        |         |
|--------|--------|---------|
| 60 KVP | 80 KVP | 100 KVP |
| 56.5%  | 47.1%  | 46.1%   |

- 30" thick textured finger/palm
- Flexibility, dexterity, and sensitivity to maintain the tactile feel
- Direct exposure attenuation rates



## Contact Us

For complete catalog, pricing and information on our custom products and kits, please contact Epimed Customer Service or your local Epimed Pain Representative. ☎ 972 373-9090. 🌐 [www.epimed.com](http://www.epimed.com)

**Protective Aprons Warranty:** Epimed International warrants that all x-ray protective aprons manufactured under its label to be free of defects in material and workmanship for a period of 2 years from the date of purchase. Seller's obligation is limited exclusively to the repair or replacement of the defective product. This warranty does not cover damage as a result of abuse, abrasion, or failure to follow care and use instructions. Customer mishandling of products is not covered under warranty. These include, but are not limited to, the following:

- Melting of apron due to heat exposure
- Puncture holes through the fabric and core materials of apron due to a sharp object
- Submersion of apron for cleaning purposes

Seller hereby disclaims all other warranties, either expressed or implied. There are no warranties which extend beyond the time frame stated. In no case shall the liability of Epimed International exceed the purchase or replacement price of the product. UltraFlex™, Secure Shield™, Secure Touch® are trademarks of Barrier Technologies®. Velcro® is a registered trademark of Velcro Industries B.V.



# Radiation Protection

## Aprons, Gloves and Accessories

Minimizing Risk from Scatter Radiation



Minimizing Risk.  
Every Procedure.  
Every Time.

# Aprons

Epimed offers a complete line of innovative radiation protection aprons to mitigate the risk of X-Ray exposure. They are made in the USA to ensure the highest standards and quality control. Our easy-to-wear aprons offer lasting radiation protection. Custom-made aprons come with a broad range of options to fit YOUR needs.



**Lightweight**



**Flexible**



**Anti-bacterial**



**Stain Resistant**

## Ready-To-Go



- Ships within 24hrs of order confirmation
- ULTRAFLEX™ anti-bacterial stain resistant fabric
- MAGNAGUARD™ magnetic thyroid collar included
- Choose from 2 styles & 4 colors



## The Flex

- Excellent weight distribution & comfort
- Two fully adjustable crossing Velcro straps
- Adjustable elastic inserts in the rear flaps
- Shoulder pads & thyroid collar included
- ULTRAFLEX™ fabric
- Choose from 11 colors




## Flex Vest and Skirt

- Elastic back for greater comfort
- Even weight distribution reduces fatigue
- Optional elastic insert in skirt
- Shoulder pads & thyroid collar included
- ULTRAFLEX™ fabric
- Choose from 11 colors


## Thyroid Collars

- 100% attenuation @ 100kvp
- Complete protection to the neck and sternum
- Disposable and visor options
- Machine washable
- Available separate, sewn, or tethered to your apron



## Protective Sleeves

- 100% attenuation @ 100kvp
- Available in full and half length
- Protects arm against scatter radiation
- Effectively shields long bones



## Personalization Options

To ensure that users find the perfect apron for their needs, different customization options are available.

### Protections Performance

Two choices of protection performance are available, related to the thickness.

**.35mm**  
offer 97% attenuation  
@ 100kvp

**.5mm**  
offer 100% attenuation  
@ 100kvp

### Core Materials

Two choices of core materials are available.

**PB FREE™**  
Lead and Toxin Free  
Super Lightweight

**PB STANDARD™**  
Traditional Lead  
Regular Weight

### Embroidery

Personalize your apron with a custom design or logo on the pocket.



### Color Choices

Apron

Ready-To-Go options



Black



Charcoal Grey



Navy



Royal Blue



Red



Pink



Yellow



Aqua



Teal



Camel



Purple

### Binding and Inside linings

Ready-To-Go



Black



Red



Silver

### Apron Sizing



Aprons are available in Standard Sizes and Custom Sizes. To facilitate your measurement, download the mobile app: **X-Ray Apron Sizer**. A helpful tool for getting precise body measurements in a few minutes. To start the measurement, use the affiliate sharing code **Epimed123**.



Your measurement  
in minutes

### Apron Tracking



The mobile app: **BarrierTrack** simplifies the management of your radiation protection garments. Track inspection histories, inventory status, and more with customized reports that can be printed or downloaded. By utilizing these reports, detailed information can be quickly and neatly organized for routine audits from authoritative bodies. **2 Year Warranty Included**



Simplify Your  
Inspection Process

## Other Accessories

In addition to the protective aprons, various accessories are offered depending on your needs.

• Wall Barriers and Shields

• Apron Racks

• Table Skirts

• Apron Cleaning Wipes

